

INITIATE AN INTRAOSSEOUS INFUSION (FAST1®)

For use of this form see TC 8-800; the proponent agency is TRADOC.

TABLE: III
REFERENCE: STP 8-68W13-SM-TG, Task: 081-833-0185, Initiate a FAST 1®

PRIVACY ACT STATEMENT

AUTHORITY: 10 U.S.C. § 3013 Secretary of the Army; AR 350-1, Army Training Leadership and Development.
PRINCIPAL PURPOSE: To ensure that accomplishment of training is properly credited to the correct individual for NREMT certification IAW AR 40-68, AR 220-1 and AR 350-1.
ROUTINE USES: Used by Unit personnel to monitor training. The DOD "Blanket Routine uses" set forth at the beginning of the Army's compilation of system of records notices may apply to this system.
DISCLOSURE: Voluntary. Failure to provide your name may result in a loss of credit for accomplishing the training or error in processing applicable favorable personnel actions. For Official Use Only.

1. Soldier (*Last Name, First Name, MI*) **2. Date** (YYYYMMDD)

SCENARIO:
 The lead Bradley infantry fighting vehicle traveling directly in front of your vehicle has been struck by an IED. After securing the immediate area, the patients are moved to safety. Following your initial assessment, you determine that you must establish peripheral intravenous access for one of the patients in order to initiate fluids. After two attempts at initiating a peripheral intravenous (IV), you are unable to establish venous access. You must establish vascular access by a sternal intraosseous (first access for shock and trauma/FAST1®) device; you have been provided the necessary medical equipment.

GRADING SHEET

TASK	COMPLETED					
	1ST		2ND		3RD	
	P	F	P	F	P	F
3. Performance Measures						
a. Assembled the necessary equipment.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
b. Spiked the IV bag and properly prepared the IV tubing.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
c. Took body substance isolation precautions.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
d. Located the suprasternal notch.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
e. Cleaned site with alcohol.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
f. Placed the target patch at landmark.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
g. Rechecked the location of the target patch.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
h. Placed the bone needle cluster into the target zone of the target patch; maintained perpendicular aspect of the introducer to the manubrium.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
i. Applied firm, increasing pressure along the axis of the introducer until a distinct release was felt and heard.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
j. Gently removed the introducer by pulling straight back.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
k. Connected the infusion tube to the right angle connector on the target patch.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
l. Connected the syringe and flushed the infusion tube with 5 ml of sterile saline.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
m. Connected the IV infusion tubing; opened IV line.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
n. Attached the protective dome to target patch and secured with tape.	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
o. Documented treatment on DD Form 1380 [Tactical Combat Casualty Care (TCCC) Card].	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

4. Demonstrated Proficiency Yes No

5. Start Time **6. Stop Time** **7. Initial Evaluator**

This form was prepared by U.S. Government employees for use in the 68W MOS. Although it contains, in part, copyrighted material from National Registry of Emergency Medical Technicians, Inc. (NREMT) skill sheets ©2011, this form has neither been prepared nor approved by NREMT. Use is restricted to guidelines contained in the Preface to TC 8-800.

INITIATE AN INTRAOSSEOUS INFUSION (FAST1®)
Grading Sheet (*cont'd*)

8. Start Time	9. Stop Time	10. Retest Evaluator
11. Start Time	12. Stop Time	13. Final Evaluator

14. Evaluator's Comments